




— :



Live time : 281 s Processing Time : P2
 XGT Dia. : 3 mm X-ray tube vol. : 50 kV
 Current : 0.460 mA
 X-ray Filter : 5 Elements Cell : Nonexistence
 Second Filter : OFF

Quant. Corr. : Standard-less

Elem.	Line	Mass [%]	2sigma [%]	Atomic [%]	Intensity [cps/mA]
24	Cr K	0.00	0.10	0.00	0.00
29	Cu K	23.35	0.10	30.29	1307.53
30	Zn K	22.47	0.09	28.34	1726.35
47	Ag K	54.16	0.12	41.37	6355.32
48	Cd K	0.00	0.02	0.00	0.34
80	Hg L	0.00	0.04	0.00	0.00
82	Pb L	0.02	0.03	0.01	0.61